

CRITTER CRADLE® CC-PRONE

STEP-BY-STEP INSTRUCTIONS



INSTRUCTIONS

Safety First!

Always practice safety for your patients and staff. Wear protective gloves, smock and eye protection. Remove jewelry before making forms. Always follow directions!

Critter Cradle[®] CC-Prone

Each kit includes a set of foaming agents, base, V-tech pieces and a special blue polyform bag.



Assemble the v-tech pieces to the base to desired width. Place the form into the bag. Mix the two bottles of foam together and vigorously shake for 15 timed seconds. Remove the cap and pour the creamy foam over the form while collapsing the accordion style bottle, leaving as little hang-up as possible.



Once foam begins to rise, carefully unfold the bag and gently press against the foam. Use your hand to spread as much air as possible out of the bag. Folding the ends inward, carefully lift the form and fold bag underneath. Spread foam evenly with your hands.



When foam is spread evenly and is rising in the mold, carefully place patient on top of the mold. Move patient around until they are centered in the form with the head and torso in a 90 degree angle from the mold. Carefully push foam in place firmly against the patient.



As foam rises around patient, it may lift the head and extremities out of the mold. Carefully press those anatomical sites deep into the foam and hold down for 15 seconds. This will set the foam and keep the anatomy from rising again. The form should be ready to use in 15 minutes.



Lift your patient out of the mold by lifting their legs and head first. Supporting the torso, lift the rest of the patient out of the mold. For larger patients, use 2 people to lift patient safely out of the mold. To place patient back in mold, lower legs and torso first, then head and front legs.



CRITTER CRADLE® BRAND CC-PRONE HELPFUL HINTS

STRENGTHENING FORMS: After pouring the foam around the edges of the form and placing the form deep into the bag, remove as much air as you can before folding the excess bag underneath the form. Allow the foam to start to rise above the form, then push down, keeping the foam level with the form, forcing it tightly against the patient

CREATING TIGHT, SECURE FORMS: Do not let the patient move. Keep constant pressure on the patient's anatomy that rises out of the foam, applying gentle yet firm pressure.

The form is useless unless you make a deep impression of the patient in the desired set-up position. Sometimes this may mean putting pressure on top of the patient's shoulders, keeping them deep in the form, allowing the foam to rise to the top of the form.

Foam that rises above the white form is useless foam. It is foam that can be diverted to another area where more foam may be needed or it can be used to make a tight fitting mold by pressing it down and holding it down for several minutes, until the form begins to cool.

CURING EACH MOLD: It is critical that the patient does not move for a full 15 minutes to allow the form to harden properly. As the form hardens, you can feel it cool to room temperature within 15 minutes from the starting point (depending upon temperature and atmospheric conditions in the room). Do not allow the patient to get up until the cradle is cool (room temperature) to the touch. If the area underneath the patient is still warm, allow it to cool and harden properly before use. This will result in a secure, comfortable tool to quickly reposition your patient each time they come in for treatment.



**CRITTER
CRADLES**
We Cradle Your Pet in Love